# LC 6000

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## Pompa SCION 6100

### Dimensioni e pesi

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No.	LED Status	Vacuum Pump State
1.	Orange – Blinking	Power On – Running
2.	Green – Blinking	Vacuum degree reached the upper limit of the control range
3.	Green – Solid	Vacuum degree reached a setting control point
4.	Orange – Solid	Error – vacuum pump has stopped

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Table 2: Mobile phase compatibility

No.	Solvent	Remarks
1.	Acid	High concentration > 1 mol/L – corrodes chamber
2.	Alkali	High concentration > 1 mol/L – corrodes chamber
3.	Salts	Salts and a high buffer solution require a wash with distilled water to reduce salt deposition as it may cause damage to the instrument
4.	Buffer Solution	Salts and a high buffer solution require a wash with distilled water to reduce salt deposition as it may cause damage to the instrument

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External I/O:

- i) RS-422
  - (a) BUSY OUT
  - (b) ERROR IN/OUT
  - (c) GRADIENT START IN
- ii) Morsetti sul retro:
  - (a) EVENT1 - 4 (Limit out)
  - (b) PUMP ON IN
  - (c) PUMP OFF IN
  - (d) PRESS OUT

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### Six-Channel Degassing Unit

Flow Path:	6 (Pump: 4 flow paths; Autosampler: 2 flow paths)
Max Flow Rate:	5.0 mL/min
Recommended Flow Rate:	<3.0 mL/min
Chamber Volume (480 µL, inlet tube not included)	
Max Pressure:	0.2 MPa* Flow path in a degassing unit should always maintain constant negative pressure.
Wetted Part Material:	Teflon® AF, PEEK, FEP, PPS (glass included).

### Instrument Status

Degasser pressure sensor error

### Auto Purge Valve Unit

Optional valve for purging automatically

Max pressure:	60MPa
Settable Flow Range:	0.500-5.0mL/min
Settable Time:	1-30 min
Wetted Materials:	Stainless Steel Stator and Steel Rotor Seal*

\*Vespe® is used in the auto purge valve. Avoid using strong acids or bases, as these may cause corrosion.

### Plunger Wash Pump

Flow Rate:	1mL/min
Delivering Solvent:	Distilled water only
Max Pressure:	65kPa

### Dynamic Mixer

Max Pressure:	40kPa
Mixer Volume:	2000µL
Mixing Method:	Stirring Bar
Recommended Flow Rate:	<2.0mL/min

### Conventional Mixer

Max Pressure:	60 MPa
Mixer Volume:	700 µL
Recommended Flow Rate:	0.4 - 1.8 mL/min

### Semi Micro Mixer

Max Pressure:	60 MPa
Mixer Volume:	200 µL
Recommended Flow Rate:	< 0.4 mL/min

#### SCION Instruments HQ

Suite 34, Livingston Business Centre  
Kirkton Road South, Livingston,  
EH54 7FA, Scotland, UK.  
Tel: +44 1506 300 200  
sales-eu@scioninstruments.com  
www.scioninstruments.com

#### SCION Instruments NL BV

Amundsenweg 22-24,  
4462 GP Goes,  
The Netherlands  
Tel: +31 (0) 113-287 600  
sales-eu@scioninstruments.com  
www.scioninstruments.com